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## Symbol Descriptions: English

	Indicates conditions or practices that could lead to illness, injury, or death ("WARNING").
	Indicates conditions or practices that could damage the equipment or other property ("CAUTION").
	Consult instructions for Use
	Do not re-use, Single use device
	Do not re-sterilize
	Catalogue Number
	Lot Number
	Use By Date
	Manufacture Date
	Quantity
	Manufacturer's contact information
	Not made with natural rubber latex
	Fragile
	Temperature limit
	Humidity limitation
	For Professional Use
	Do not use if package is damaged
	Sterilized using ethylene oxide (only applicable to versions marked as sterile on unit packaging)
	TYPE BF APPLIED PART

### Indications for Use

The single-use lighted retractator is indicated for use to retract soft tissue in order to enhance exposure during surgical procedures.

### ⚠️ Warnings and Cautions

**WARNING:** Inspect each device for damage prior to use  
**WARNING:** When in transit or storage, device may be subject to damage beyond the control of the manufacturer or supplier.  
**WARNING:** Never use device with laser equipment.  
**WARNING:** Do not allow conductive instruments to contact LED's.  
**WARNING:** Treat used devices as bio hazardous infectious material. Dispose of used device in suitable disposal unit or in accordance with local regulations.  
**WARNING:** Devices are not compatible with gamma radiation or autoclave sterilization and should never be re-sterilized by any method.  
**WARNING:** No modification of this equipment is allowed  
**WARNING - Risk of reuse:**  
The strength of the device may be significantly compromised and the resulting risk to the patient and/or user is unknown.  
The light may not function properly.  
**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician.  
**CAUTION:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ONETRAC, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.  
**CAUTION:** Electromagnetic disturbances may create imperceptible fluctuations in the intensity of the LED light.


### Warranty

Please visit [www.obpmedical.com](http://www.obpmedical.com) "Terms and Conditions" page for warranty information.  
OBP Medical should be notified immediately via phone or email of any defective units.


**REF** C100110, C100140, C100150  
C100165, C100170, C100180  
C100190, C100300

**OBP Medical**  
360 Merrimack Street  
Building 9  
Lawrence, MA 01843 USA  
+1-978-291-6853  
[www.obpmedical.com](http://www.obpmedical.com)

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