

ONETRAC™

REF C100110, C100140, C100150, C100165, C100170, C100180, C100190

English

- Do not reuse, single-use device
Catalog number
Lot number
Use by date
Date of manufacture and country of manufacture
Quantity
Manufacturer
Fragile during transport and storage
Temperature storage limit
Humidity storage limit
Do not use if package is damaged and consult IFU
Single-sterile barrier system with protective packaging inside, sterilized using ethylene oxide
Type BF applied part, internally powered ME equipment
Authorized representative in the European Community/ European Union
Keep dry during transport and storage
Medical device
Unique device identifier
Size/model
Consult Instructions for Use
Do not resterilize
Vertically falling drops; no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical

Product overview:

ONETRAC is a single-use cordless retractor with integrated LED light source. ONETRAC is a handheld instrument that provides soft tissue retraction during surgery (essential performance), provides illumination to the surgical site when the integrated lighting feature is turned on, and provides smoke evacuation from the surgical field when attached to a vacuum source.

Intended Use/Purpose and Indications for Use:

The single-use lighted retractor is indicated for use to retract soft tissue in order to enhance exposure during surgical procedures.

Instructions Prior to Use

- Verify that the EtO indicator sticker on the unit box is BLUE - do not use if EtO indicator sticker on unit box is RED.
Inspect sterile barrier for damage - if the sterile barrier packaging is damaged do not use and source a new device.
Inspect each device for damage - if device is damaged do not use and source a new device.
Confirm expiry date prior to use - if current date is beyond the expiry date, do not use and source a new device.
Confirm device functionality according to instructions for use - if device is defective, do not use and source a new device.

obp Medical
360 Merrimack St., Bld. 9
Lawrence, MA 01843 USA
www.obpsurgical.com
sales@obpsurgical.com
+1 (978) 291-6853

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WARNING: If EtO indicator sticker is RED or if the sterile barrier packaging is damaged, contents may not be sterile. Do not resterilize.
WARNING: If used beyond expiry date, contents may no longer be sterile. Do not resterilize.
WARNING: Use of a damaged device may result in patient or user injury.

Instructions For Use

- 1. Lighting: Before using the device, confirm LED(s) on the blade turn on by sliding the black on/off switch (A) forward to the “on” position.
2. Smoke Evacuation (available on all blade sizes except 55x8): If smoke evacuation is desired attach smoke evacuation tubing (not included) by press-fitting the female tubing connector onto the smoke evacuation port (B) on the handle.
3. Optional Battery Disposal After Use: If alkaline button batteries require separate disposal from the device, locate the small battery cover on the end cap next to the smoke evacuation port.
4. Device Disposal After Use: Treat used device as biohazardous infectious material and dispose of in accordance with state and local regulations and/or facility procedures.

WARNING: Exercise caution while dissecting around the device. Direct contact with an electro-surgical instrument may cause physical damage to the device or the patient.
WARNING: Applying excessive force to tissue with APPLIED PART (D) may result in tissue damage and/or device breakage.
WARNING: Unauthorized modification of the device (for example - re-processing, replacing battery, or re-sterilizing) may cause damage to the battery resulting in chemical burns or fire/explosion and may result in a non-sterile device.
WARNING: Improper disposal of used device may lead to spread of infection.
WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.
CAUTION: Do not submerge device in fluid, as this may result in loss of light function.
CAUTION: In case of loss of illumination during use, retrieve alternative light source.
CAUTION: Surgical smoke may pose a health hazard.
Conforms to IEC 60601-1:2005, IEC 60601-1-2:2014, INTERNALLY POWERED ME EQUIPMENT; CISPR 11 CLASS A, GROUP 1

Warranty Disclaimer

obp products are provided “as is” and without warranties of any kind, whether express or implied. To the fullest extent permissible pursuant to applicable law, obp disclaims all warranties, express or implied, including, but not limited to, implied warranties of merchantability and fitness for a particular purpose and non-infringement.

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