# ONETRAC™ LX

C100210, C100220, C100230, C100240, C100260, C100270, C100280, C100290

## English

Contains two (2) non-rechargeable Li-metal batteries. Dispose of batteries in accordance with local hazardous waste regulations/requirements.

Do not reuse, single-use device

Catalog number

Lot number

Use by date

Date of manufacture and country of manufacture

Quantity

Manufacturer

Fragile during transport and storage

Temperature storage limit

Humidity storage limit

Do not use if package is damaged and consult IFU

Sterilized using ethylene oxide

Type BF applied part, internally powered ME equipment

Keep dry during transport and storage

Medical device

Unique device identifier

Size/model

Consult Instructions for Use

Do not resterilize

Single sterile barrier system

Single sterile barrier system with protective packaging outside

Authorized representative in the European Community/ European Union

Vertically falling drops; no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical

IPX2

### Product overview:

ONETRAC LX is a single-use cordless retractor with integrated multi-LED light source. ONETRAC LX is a handheld instrument that provides soft tissue retraction during surgery (essential performance), provides illumination to the surgical site when the integrated lighting feature is turned on, and provides smoke evacuation from the surgical field when attached to a vacuum source. ONETRAC LX is available in several blade sizes. It is a single-use device provided sterile and is not intentionally manufactured using natural rubber latex. Once the device is illuminated, the device provides a minimum of 180 minutes of continuous illumination. The device is intended for use only by qualified medical personnel possessing professional training in the surgical procedures being performed, while in a surgical suite/operating room with a maximum operating temperature (local ambient temperature) up to 25°C. The indirect clinical benefit intended is enhanced visibility in the surgical pocket.

### Intended Use/Purpose and Indications for Use:

The single-use lighted retractor is indicated for use to retract soft tissue in order to enhance exposure during surgical procedures.

### Instructions Prior to Use

- 1. Verify that the EIO indicator sticker on the unit box is BLUE - do not use if EIO indicator sticker on unit box is RED.
2. Inspect sterile barrier for damage - if the sterile barrier packaging is damaged do not use and source a new device.
3. Confirm activator tab (A, F) is only partially inserted into the handle while device is sealed in the sterile barrier - if activator tab is depressed into the handle within the sealed sterile packaging do not use and source a new device.
4. Inspect each device for damage - if device is damaged do not use and source a new device.
5. Confirm expiry date prior to use - if current date is beyond the expiry date, do not use and source a new device.
6. Confirm device functionality according to instructions for use - if device is defective, do not use and source a new device.

WARNING: If EIO indicator sticker is RED or if the sterile barrier packaging is damaged, contents may not be sterile. Do not resterilize.
WARNING: If used beyond expiry date, contents may no longer be sterile. Do not resterilize.
WARNING: Use of a damaged device may result in patient or user injury.

## Instructions For Use

### Single-blade Retractor Instructions

- 1. Battery Activation: Depress colored battery activator tab (A) into handle, do not remove from handle.
2. Lighting: Before using the device, confirm all LEDs on the blade turn on by pressing the black on/off button (B) 3x.
3. Smoke Evacuation: If smoke evacuation is desired, attach smoke evacuation tubing (not included) by press-fitting the female tubing connector onto the smoke evacuation port (C) on the handle.
4. Battery Disposal After Use: Prior to disposal of device, with battery door facing downward, open the battery door by pulling on the latch (D).
5. Device Disposal After Use: Once batteries are removed, treat used device as biohazardous infectious material and dispose of in accordance with state and local regulations and/or facility procedures.

### Dual-blade Retractor Instructions

- 1. Battery Activation: Depress colored battery activator tab (F) into handle, do not remove from handle.
2. Lighting: Before using the device, confirm all LEDs on both blades turn on: Press the black rocker button (G) 2x closest to the blades you are activating.
3. Smoke Evacuation: If smoke evacuation is desired: Insert narrower side of included smoke evacuation tubing adapter into the handle port (I) on the opposing side to the blade being used.
4. Battery Disposal After Use: Prior to disposal of device, with battery door facing downward, open the battery door by pulling on the latch (J).
5. Device Disposal After Use: Once batteries are removed, treat used device as biohazardous infectious material and dispose of in accordance with state and local regulations and/or facility procedures.

WARNING: Use of alternative smoke evacuation tubing adapters may result in loss or reduced suction and therefore reduced visibility in the surgical cavity during cautery and exposure to fumes from surgical smoke.

WARNING: Exercise caution while dissecting around the device. Direct contact with an electrosurgical instrument may cause physical damage to the device or the patient.

WARNING: Applying excessive force to tissue with APPLIED PART (E, K) may result in tissue damage and/or device breakage. In case of device breakage during use, ensure no fragments or batteries remain in patient and retrieve a new device.

WARNING: Unauthorized modification of the device (for example-reprocessing, replacing battery or re-sterilizing) may cause damage to the battery resulting in chemical burns or fire/explosion and may result in a non-sterile device.

WARNING: Improper disposal of used device may lead to spread of infection.
WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ONETRAC LX. Otherwise, degradation of the performance of the ONETRAC LX could result.

WARNING: Do not remove sterile device tray from protective unit box until use. Preemptive removal may result in damage to sterile barrier which may cause patient infection.

CAUTION: Do not submerge device in fluid, as this may result in loss of light function.

CAUTION: In case of loss of illumination during use, retrieve alternative light source.

CAUTION: Possibly hazardous optical radiation emitted from the LEDs. Do not stare at the LEDs while on as this may be harmful to the eye.

CAUTION: In extreme cases, electromagnetic disturbances may cause the light to turn OFF. To turn the light back ON, press the ON/OFF button.

CAUTION: Surgical smoke may pose a health hazard. In case of loss of smoke evacuation during use, detach smoke evacuation tubing from device and use tubing directly.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Conforms to IEC 60601-1:2005, IEC 60601-1-2:2014, INTERNALLY POWERED ME EQUIPMENT; CISPR 11 CLASS A, GROUP 1

### Warranty Disclaimer

obp products are provided "as is" and without warranties of any kind, whether express or implied. To the fullest extent permissible pursuant to applicable law, obp disclaims all warranties, express or implied, including, but not limited to, implied warranties of merchantability and fitness for a particular purpose and non-infringement. Notwithstanding the aforementioned, obp will warranty against manufacturer defects for single-use devices used prior to their expiration date. The above limitations are in full effect within states where such warranties are permitted.

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