
# ONETRAC™ LXS Yankauer

**REF C110200, C110210, C110220**

## English

- Contains two (2) non-rechargeable Li-metal batteries. Dispose of batteries in accordance with local hazardous waste regulations/requirements.
- Do not reuse, single-use device
- Catalog number
- Lot number
- Use by date
- Date of manufacture and country of manufacture
- Quantity
- Manufacturer
- Fragile during transport and storage
- Temperature storage limit
- Humidity storage limit
- Do not use if package is damaged and consult IFU
- Sterilized using ethylene oxide
- Type BF applied part, internally powered ME equipment
- Keep dry during transport and storage
- Medical device
- Unique device identifier
- Size/model
- Consult Instructions for Use
- Do not resterilize
- Single sterile barrier system
- Single sterile barrier system with protective packaging outside
- Authorized representative in the European Community/ European Union
- IPX2

**Product overview:** ONETRAC LXS Yankauer is a single-use suction device with integrated cordless radial LED light source. ONETRAC LXS Yankauer is a handheld instrument that provides aspiration from the surgical field (essential performance) when attached to a vacuum source and provides illumination to the surgical site when the integrated lighting feature is turned on. The ONETRAC LXS Yankauer is available in multiple tip configurations and is provided with a stylet. It is a single-use device provided sterile and is not intentionally manufactured using natural rubber latex. Once the device is illuminated, the device provides a minimum of 180 minutes of continuous illumination. The device is intended for use only by qualified medical personnel possessing professional training in the surgical procedures being performed, while in a surgical suite/operating room with a maximum operating temperature (local ambient temperature) up to 25° C. The indirect clinical benefit intended is enhanced visibility in the surgical pocket.

**Intended Use/Purpose and indications for Use:** Single-use, disposable hand tool with an integrated LED light used for sub cutaneous suction in surgical procedures.

- Instructions Prior to Use**
1. Inspect sterile barrier for damage – if the sterile barrier packaging is damaged do not use and source a new device.
  2. Confirm activator tab (A) is only partially inserted into the handle while device is sealed in the sterile barrier - if activator tab is depressed into the handle within the sealed sterile packaging do not use and source a new device.
  3. Confirm stylet is packaged inside the sterile barrier with the device - if stylet is not packaged with the device, do not use and source a new device.
  4. Inspect each device for damage - if device is damaged do not use and source a new device.
  5. Confirm expiry date prior to use - if current date is beyond the expiry date, do not use and source a new device.
  6. Confirm device functionality according to instructions for use - if device is defective, do not use and source a new device.

WARNING: If EtO indicator sticker is RED or if the sterile barrier packaging is damaged, contents may not be sterile. Do not resterilize.  
WARNING: If used beyond expiry date, contents may no longer be sterile. Do not resterilize.  
WARNING: Use of a damaged device may result in patient or user injury.

## Deutch

- Instructions For Use**
1. **Battery Activation:** Depress colored activator tab (A) into handle, do not remove from handle.
  2. **Lighting:** Before using the device, confirm LED turns on by pressing the black on/off button (B). If illumination is intermittent or does not turn on, retrieve a new device. When device illumination slowly turns off and on autonomously, the batteries are fully depleted. Turn off light when not in use to preserve battery life.
  3. **Suction:** Attach suction tubing (not included) by press-fitting the female tubing connector onto the suction port (C) on the handle. Confirm the tubing is securely attached to the port (C), and confirm vacuum power at distal suction tip before use. Use included stylet to unclog rigid suction cannula if necessary. **WARNING:** If device cannot be unclogged, source a new device.
  4. **Battery Disposal after use:** Prior to disposal of device, with battery door facing downward, open the battery door by pressing on the latch (D) and by lifting the battery door away from the handle. Continue to lift the battery door until the batteries are sitting within the battery door. Tilt the battery door until the batteries fall out of the battery door and into the battery disposal container. Dispose of lithium primary batteries according to state and local regulations and/or facility procedures.
- WARNING: Exercise caution while dissecting around the device. Direct contact with an electrosurgical instrument may cause physical damage to the device or the patient.  
WARNING: Applying excessive force to tissue with APPLIED PART (F) may result in tissue damage and/or device breakage. In case of device breakage during use, ensure no fragments or batteries remain in patient and retrieve a new device.  
WARNING: Unauthorized modification of the device (for example-reprocessing, replacing battery, or re-sterilizing) may cause damage to the battery resulting in chemical burns or fire/explosion and may result in a non-sterile device.  
WARNING: Improper disposal of used device may lead to spread of infection.  
WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.  
WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ONETRAC LXS. Otherwise, degradation of the performance of the ONETRAC LXS could result.  
WARNING: Do not remove sterile device tray from protective unit box until use. Preemptive removal may result in damage to sterile barrier which may cause patient infection.  
CAUTION: Do not submerge device in fluid, as this may result in loss of light function.  
CAUTION: In case of loss of illumination during use, retrieve alternative light source.  
CAUTION: Possibly hazardous optical radiation emitted from the LEDs. Do not stare at the LEDs while on as this may be harmful to the eye.  
CAUTION: In extreme cases, electromagnetic disturbances may cause the light to turn OFF. To turn the light back ON, press the ON/OFF button.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
- Conforms to IEC 60601-1:2005,IEC 60601-1-2:2014, INTERNALLY POWERED ME EQUIPMENT; CISPR 11 CLASS A, GROUP 1
- Warranty Disclaimer**
- obp products are provided "as is" and without warranties of any kind, whether express or implied. To the fullest extent permissible pursuant to applicable law, obp disclaims all warranties, express or implied, including, but not limited to, implied warranties of merchantability and fitness for a particular purpose and non-infringement. Notwithstanding the aforementioned, obp will warranty against manufacturer defects for single-use devices used prior to their expiration date. The above limitations are in full effect within states where such warranties are permitted.

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