


ONETRAC[™] LX D

Single-Use Cordless Dual-Blade Retractor with Integrated Multi-LED Light Source



**CORDLESS ILLUMINATION
WITH INTEGRATED SMOKE
EVACUATION**

ONETRAC LX D is a single-use, cordless dual-blade surgical retractor with integrated multi-LED light source and dual smoke evacuation channels.

This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the surgical setting by making reusable devices and fiber-optic light cables obsolete.

Illumination and smoke evacuation in the depths of a surgical pocket can be met within seconds – no assembly, additional parts/components, or fiber-optic cables required.

To learn more visit www.obpsurgical.com/lx

TO ORDER

EMAIL: SALES@OBPSURGICAL.COM

CALL: 1-978-291-6853


WEB: OBPSURGICAL.COM

- Item # C100270 (Dual-Blade)
 - 50mm x 20mm
 - 30mm x 15mm
- Sold individually (1 EA/BX)
- Each unit comes packed sterile

ABOUT US

obp is the leading global developer of single-use, self-contained, illuminating surgical devices. We integrate light with device — making them one — creating effective and simplified single-use solutions. Our innovative products are used in surgery centers and hospitals throughout the U.S., as well as healthcare facilities worldwide.

CUT THE CORD



PRODUCT FEATURES

- Disposable, single-use device reduces risk of cross-contamination
- Integrated cordless LED lights are cooler, providing direct illumination and peak power for over 3 hours of continuous use
- Patent-pending dual-blade LED circuitry enables instantaneous transition between LEDs on either blade with the press of a button
- Patent-pending dual smoke evacuation channels allow for efficient mid-procedure transition from one blade to another (custom smoke evacuation adapter included with each unit)
- Color-coded battery activation tab allows user to easily identify and differentiate blade sizes
- Newly designed ergonomic handle for improved maneuverability
- Atraumatic grip at distal end of longer blade to maximize tissue adherence
- Patented touch-free battery removal allows for disposal of batteries without battery contact for easy recycling

GPO CONTRACTS

PREMIER CONTRACT PP-OR-1741

HEALTHTRUST CONTRACT 39624



360 Merrimack Street, Lawrence, MA 01843 USA
+1.978.291.6853 | www.obpsurgical.com

PL-10-0915 D

