

# ONETRAC™ LXS

Single-Use Suction Device with Integrated Cordless Radial LED Light Source



SINGLE-USE CORDLESS  
ILLUMINATION

**ONETRAC LXS is a single-use suction device with integrated, cordless radial LED light source.**

This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the surgical setting by making reusable devices and fiber-optic light cables obsolete.

**Illumination and aspiration** in the depths of a surgical pocket can be met within seconds – no assembly, additional parts/components, or fiber-optic cables required.

To learn more visit [www.obpsurgical.com/lxs](http://www.obpsurgical.com/lxs)

## TO ORDER

EMAIL: [SALES@OBPSURGICAL.COM](mailto:SALES@OBPSURGICAL.COM)

CALL: 1-978-291-6853


WEB: [OBPSURGICAL.COM](http://OBPSURGICAL.COM)

- Item # C110200 (Crown Tip)
- Sold individually (1 EA/BX)
- Each unit comes packed sterile

## ABOUT US

**obp** is the leading global developer of single-use, self-contained, illuminating surgical devices. We integrate light with device — making them one — creating effective and simplified single-use solutions. Our innovative products are used in surgery centers and hospitals throughout the U.S., as well as healthcare facilities worldwide.

# CUT THE CORD



## PRODUCT FEATURES

- Disposable, single-use device reduces risk of cross-contamination
- Integrated cordless LED lights are cooler, providing direct illumination and peak power for over 3 hours of continuous use
- Radial light coverage allows for full visibility of a natural cavity or surgical pocket
- Patent pending LED driven lens & reflector technology enable an ultra-bright spotlight housed in a low-profile package
- Centrally located on/off button
- 4.5mm cannula channel allows for rapid suction
- Color-coded battery activation tab
- Curved cannula for easy access and targeted placement within a deep pocket
- Ergonomic handle for improved maneuverability
- Atraumatic clog resistant crown tip at distal end
- Patented touch-free battery removal allows for disposal of batteries without battery contact for easy recycling

## GPO CONTRACTS

PREMIER CONTRACT PP-OR-1741

HEALTHTRUST CONTRACT 39624



360 Merrimack Street, Lawrence, MA 01843 USA  
+1.978.291.6853 | [www.obpsurgical.com](http://www.obpsurgical.com)

PL-11-0901 D

