


# ONETRAC™ LXS

Single-Use Suction Device with Integrated Cordless Radial LED Light Source



## Frazier

**ONETRAC LXS Frazier is a single-use suction device with integrated, cordless radial LED light source.**

This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the surgical setting by making reusable frazier suction devices and associated fiber-optic light cables obsolete.

Illumination in the depths of a surgical pocket can be met within seconds – no assembly, additional parts/components, or fiber-optic cables required.

## TO ORDER

---

EMAIL: [SALES@OBPSURGICAL.COM](mailto:SALES@OBPSURGICAL.COM)

CALL: 1-978-291-6853

WEB: [OBPSURGICAL.COM](http://OBPSURGICAL.COM)

- Item # C110310 (8 Fr)
- Item # C110320 (12 Fr)
- Sold individually (1 EA/BX)
- Each unit comes packed sterile
- Stylet included with each unit

## ABOUT US

obp is the leading global developer of single-use, self-contained, illuminating surgical devices. We integrate light with device — making them one — creating effective and simplified single-use solutions. Our innovative products are used in surgery centers and hospitals throughout the U.S., as well as healthcare facilities worldwide.

# CUT THE CORD




## ALSO AVAILABLE

**ONETRAC LX S Yankauer** single-use suction device with integrated, cordless radial LED light source

Item # C110200 (Crown Tip)

Item # C110220 (Taper Tip)

Item # C110210 (Metal Tip)



## PRODUCT FEATURES

- Disposable, single-use device reduces risk of cross-contamination
- Integrated cordless LED lights provide direct illumination for over 3 hours of continuous use
- Teardrop style flow control
- Centrally located on/off button is recessed to prevent accidental contact and preserve battery life when not in use
- Color-coded battery activation tab
- Curved cannula for easy access and targeted placement within a narrow pocket
- Newly designed ergonomic handle
- Patented touch-free battery removal allows for disposal of batteries without battery contact for easy recycling
- Available in 8 and 12 French

## GPO CONTRACTS

PREMIER CONTRACT PP-OR-2239

HEALTHTRUST CONTRACT 39624



360 Merrimack Street, Lawrence, MA 01843 USA  
+1.978.291.6853 | [www.obpsurgical.com](http://www.obpsurgical.com)

PL-11-0906 REV-D

