

# ONETRAC<sup>TM</sup> LX

Single-Use Cordless Retractor with Integrated Multi-LED Light Source



**CORDLESS ILLUMINATION  
WITH INTEGRATED SMOKE  
EVACUATION**

**ONETRAC LX is a single-use, cordless surgical retractor with integrated multi-LED light source and dual smoke evacuation channels.**

This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the surgical setting by making reusable, handheld surgical retractors and associated fiber-optic light cables obsolete.

Illumination in the depths of a surgical pocket can be met within seconds – no assembly, additional parts/components, or fiber-optic cables required.

## TO ORDER

EMAIL: [SALES@OBPSURGICAL.COM](mailto:SALES@OBPSURGICAL.COM)

CALL: 1-978-291-6853

WEB: [OBPSURGICAL.COM](http://OBPSURGICAL.COM)

- Item # C100220 (175mm x 30mm)
- Item # C100260 (135mm x 60mm)
- Item # C100210 (135mm x 30mm)
- Item # C100290 (135mm x 20mm)
- Item # C100230 (90mm x 60mm)
- Item # C100280 (90mm x 30mm)
- Item # C100240 (90mm x 20mm)
- Item # C100270 (LXD)
- Sold individually (1 EA/BX)

## ABOUT US

**obp** is the leading global developer of single-use, self-contained, illuminating surgical devices. We integrate light with device — making them one — creating effective and simplified single-use solutions. Our innovative products are used in surgery centers and hospitals throughout the U.S., as well as healthcare facilities worldwide.

# CUT THE CORD



## AVAILABLE IN 7 SIZES



Multiple LEDs for maximum illumination strength and coverage of a deep pocket

Newly designed ergonomic handle

## ALSO AVAILABLE

**ONETRAC LXD** single-use cordless dual-blade retractor with integrated multi-LED light source and dual smoke evacuation channels

Item # C100270 (LXD)



360 Merrimack Street, Lawrence, MA 01843 USA  
+1.978.291.6853 | [www.obpsurgical.com](http://www.obpsurgical.com)

PL-10-0914-J

## PRODUCT FEATURES

- Disposable, single-use device reduces risk of cross-contamination
- Integrated cordless LED lights provide direct illumination and over 3 hours of continuous use
- On/off button allows user to cycle through multiple lighting configurations and preserve battery life when not in use
- Integrated dual smoke evacuation channels for rapid smoke suction
- Color-coded battery activation tab allows user to easily identify and differentiate blade sizes
- Graded centimeter reference markings along the blade (blade size-specific)
- Newly designed ergonomic handle
- Custom blade grip at the distal tip is an alternative to traditional teeth
- Patented touch-free battery removal allows for disposal of batteries without battery contact for easy recycling

## GPO CONTRACTS

PREMIER CONTRACT PP-OR-2239

HEALTHTRUST CONTRACT 39624

